

RNA probe synthesis with Ambion Megascript SP6 or T7 Kit (30-5-08)

Set up a 20 μ l reaction with template: 2 μ g for linearized plasmid, 800ng for PCR product. Keep buffer and reaction at RT during pipetting the reaction.

	total volume	20 μ l
Megascript buffer (10x)		2
ATP (different for SP6 and T7)		2
CTP (different for SP6 and T7)		2
GTP (different for SP6 and T7)		2
UTP (different for SP6 and T7)		1.5
DIG-11-UTP (10nmol/ μ l Boehringer)		4
H ₂ O (MilliQ is fine)		variable
DNA template (max. 4.5 μ l)		variable
RNA polymerase (T7 or SP6 Megascript)		2

6 hrs @ 37 °C.

Prepare hyb buffer, RNA gel

DNase treatment

Treat with 1 μ l DNase (Megascript) for 15min at 37°C

RNA precipitation

Add 70 μ l H₂O (RNase free)

Add 10 μ l Na-Acetate 3M pH 5.5 (final 0.3M)

Add 300 μ l of ice-cold EtOH (100 %) (3 volumes)

Precipitate for 30' @ -70 °C, spin 30' @ 4°C, 15000 rpm.

Wash pellet with EtOH (70 %, DEPC H₂O), remove residual liquid and air dry (5').

Dissolve pellet in 20 μ l DEPC H₂O

Measure concentration

Take out 0.5 μ l of the sample and add 3.5 μ l RNase free H₂O (dilution 1:8)

Spec on nanodrop and run 300ng on gel

(The average total yield of the 20 μ l reaction is between 30-40 μ g)

Store 10 μ l of RNA probe at -80 °C

Fill up 9 μ l with 150 μ l Hyb buffer (depending on concentration, final conc. 50-100ng/ μ l) and store at -20 °C.